

Shotgun Lipidomics of prostate cancer cells using ESI-MS Shimadzu 8050 and simplified data analysis by SimLipid software

Rajesh GUPTA¹, Martin SADOWSKI², Ningombam Sanjib Meitei^{3,4} and Stephen J. BLANKSBY¹



¹Central Analytical Research Facility, Queensland University of Technology, Australia
²Institute for Health and Biomedical Innovation, Queensland University of Technology, Australia
³PREMIER Biosoft, Indore, India, ⁴PREMIER Biosoft, CA, USA

Introduction

Lipidomics is a relatively recent omics field of research which includes complex lipidome analysis. It is an emerging field in biomedical research as lipids play an important role in cell, tissue and organ physiology and have potential as biomarkers of disease or treatment success. Shotgun lipidomics¹ involves identification and quantification of lipids by direct infusion of complex lipid samples into the mass spectrometer without any chromatographic separation. In 2017 prostate cancer is the most commonly diagnosed cancer in Australia among males. It is estimated that the risk of a male being diagnosed with prostate cancer by his 85th birthday will be 1 in 7.

Method


LNCaP cells with different treatments




Lipid class	Ion mode	Scan type	S/N > 3
Head group scan			
PC/SM	+Ve	PI 184.1	✓
PE	+Ve	NL 141	✓
PS	+Ve	NL 185	✓
PS	-Ve	NL 87	✓
PG	+Ve	NL 172	✓
PA	+Ve	NL 98	X
PI	-Ve	PI 241	✓
Ceramides	+Ve	PI 364.3	X
Cholesterol	+Ve	PI 369.4	X
Esters			
Triacylglycerides			
TAG (14:0)	+Ve	NL 245	✓
TAG (16:1)	+Ve	NL 271	✓
TAG (16:0)	+Ve	NL 273	✓
TAG (18:3)	+Ve	NL 295	✓
TAG (18:2)	+Ve	NL 297	✓
TAG (18:1)	+Ve	NL 299	✓
TAG (18:0)	+Ve	NL 301	✓

Lipid class	Ion mode	Scan type	S/N > 3
Triacylglycerides contd.			
TAG (20:5)	+Ve	NL319	✓
TAG (20:4)	+Ve	NL 321	✓
TAG (20:3)	+Ve	NL 323	✓
TAG (20:2)	+Ve	NL325	✓
TAG (20:1)	+Ve	NL 327	✓
TAG (20:0)	+Ve	NL 329	✓
TAG (22:6)	+Ve	NL 345	✓
TAG (22:5)	+Ve	NL 347	✓
TAG (22:4)	+Ve	NL 349	✓
TAG (22:3)	+Ve	NL 351	✓
TAG (22:2)	+Ve	NL 353	✓
TAG (22:1)	+Ve	NL 355	✓
TAG (22:0)	+Ve	NL 357	✓
TAG (18:3)	+Ve	NL 295	✓
TAG (18:2)	+Ve	NL 297	✓
TAG (18:1)	+Ve	NL 299	✓
TAG (18:0)	+Ve	NL 301	✓

Example Phospholipid scans





Example Triacylglyceride scans




Good quality mass spectra obtained from different classes of lipids

High throughput data analysis

Import raw data from native file format (.lcd) into SimLipid software in batch mode



MetaboAnalyst software



Export SimLipid results into CSV file, and perform Statistical analysis using MetaboAnalyst software

Summary

High throughput data analysis was performed by directly importing raw Shimadzu data files (.lcd format) into SimLipid. Preliminary data analysis with Simlipid has shown some promising functionalities for both qualitative and quantitative analysis. The statistical analysis revealed very good correlation by PCA analysis and good stratification of different treatments.

Future work will include addition of class specific internal standards for absolute quantitation.

References:

- Han X and Gross R. W. The journal of Lipid Research, (2003) 44: 1071-1079.
- Matyash V *et al.* The journal of Lipid Research, (2008) 49: 1137-46.